



Bridge You and Nano

Exponential Business and Technologies Company

Quantitative Composition Determination of Powder Mixtures Using XRD

Material composition is an essential piece of information for engineering, quality assurance, trace element analysis, process controls, and research and development. In ceramics, alloys, steels, geology, and many other fields knowing the composition and crystallographic phase of a material is vital to accurately predicting its performance. With x-ray diffraction (XRD), it is possible to determine the type of material and weight fraction of multiple components in the same sample. XRD has the advantage of being non-destructive and does not require large sample volumes. This method works best for any polycrystalline sample with randomly oriented grains that are less than 10 μm in size, which includes many metals, alloys, geological samples, powders, ceramics, and cements.




Figure 1. Powder mixture in low-background sample holder during XRD experiment.

There are three stages to measuring composition with XRD: acquiring the diffraction pattern, identifying the crystallographic phases present in the sample, and refining a model to determine the amount of each phase. Materials are identified by comparing the sample's diffraction pattern to those in a pattern library containing hundreds of thousands of standards. Each crystallographic phase has a unique diffraction fingerprint, allowing for the identification of the materials present in the sample. The final step is fitting the whole powder pattern to a theoretical model. The total intensity of each phase's pattern is proportional to the amount present in the sample. A model

of the phases is used to calculate a theoretical diffraction pattern which includes strain, changes in stoichiometry, texture, or sample and instrument misalignments. These variables are refined using the Rietveld method and produce composition measurements good to 5 wt % for many applications.

In this demonstrative experiment, a mixture of two known, polycrystalline powders was tested for phase identification and composition measurement using a Rigaku SmartLab X-ray Diffractometer. The powder mixture was placed in a low-background sample holder and pressed to produce a smooth, flat surface, as shown in Figure 1. The Cu x-ray source, detector, and sample were all aligned. A nickel foil was used to selectively absorb the Cu $K\beta$. The powder diffraction pattern was measured with the D/TeX 1D detector at $2^\circ/\text{min}$ over a 2θ range of 15 to 110° .




Figure 2. Experimental and calculated XRD patterns of the sample. The peak positions of anatase, corundum, and rutile are marked with purple, green, and gold triangles, respectively.

The experimental diffraction pattern from the powder mixture is shown above in Figure 2. The phases were identified as corundum (Al_2O_3) and anatase (TiO_2) with the PDXL 2 analysis program using the Crystallography Open Database powder diffraction library. Rutile (another form of TiO_2) may be present at or below 1.1 wt %. Reference markers below the data indicate the ideal peak positions for each phase. The whole pattern was fit to a model in PDXL 2 in which the Lorentz polarization, the $K_{\alpha 1}$ and $K_{\alpha 2}$ splitting, and sample self-absorption were accounted for. The weight fraction, lattice constants, thermal vibration parameters variables were fit to the experimental pattern for each phase. The resulting calculation from the best-fit model closely matches the experimental pattern ($R_{\text{wp}} = 2.55\%$) and is shown as the red line in Figure 2.

As an independent check for the XRD results, the elemental composition of the sample was determined using Energy Dispersive X-ray Spectroscopy (EDS). The weight percent of each phase was calculated assuming ideal stoichiometry. As EDS does not distinguish between crystallographic phase, the weight percent of TiO_2 is a sum of anatase and rutile. As seen below in Table 1, the composition measurements from XRD and EDS are in agreement.

Table 1 Powder Mixture Composition Results from XRD and EDS

Chemical Formula	Phase	XRD (wt. %)	EDS (wt %)
Al_2O_3	Corundum	59.8 ± 2.2	62.48 ± 3.55
TiO_2	Anatase	39.1 ± 2.2	37.52 ± 4.09
	Rutile	1.1 ± 2.1	